



**Virtual Biomedical and STEM/STEAM Education**  
**2021-1-HU01-KA220-HED-000032251**

**Laboratory hands on skill course**  
**Laparoscopy VR Trainer**  
**Controller Adapter Usage Manual**



Funded by the European Union. Views and opinions expressed are however those of the author(s) only and do not necessarily reflect those of the European Union or the European Education and Culture Executive Agency (EACEA). Neither the European Union nor EACEA can be held responsible for them.



**PÉCSI TUDOMÁNYEGYETEM**  
**UNIVERSITY OF PÉCS**

**U. PORTO**



1. 3D print the file out of a softer type of PLA, no need for harder or tougher material, this is meant to be cheap and accessible.  
**Filename:** Q2controller\_laparoscopyAdapter.stl
2. Get 2 rods of preferably aluminium, plastic or wood with the following dimensions. Radius: 3 mm (Diameter: 6mm), Length: 360 mm (36 cm). We used vitrage-curtain holders and modeller plastic rods, but anything similar works.
3. (Optional) Get 2 pieces of screws with a radius of 2mm (diameter: 4mm) and attach the rod to the base of the adapter. Glue may be used as well.
4. Snap the adapter into the ring of the Quest 2 controller. PLA is flexible enough to not break. The ring of the controller is flexible enough to not break.
5. Get a cardboard box and make 2 holes in the middle, 160mm (16cm) apart from each other.

TIP: When starting the app, rest the controllers in the holes while the startup logos play.



PÉCSI TUDOMÁNYEGYETEM  
UNIVERSITY OF PÉCS

U. PORTO

